

Specification Sheet

Product Name:

Recombinant c-Myc Antibody

SKU #:

Ab097000

Clone ID:

EP121

Isotype:

Rabbit IgG

Source:

293 supernatant

Formulation:

Liquid

Concentration:

0.38 mg/mL

Buffer Solution:

10mM PBST, 0.015% Sodium azide, pH 7.4

Storage Conditions:

Store at -20°C, Avoid repeated freezing and thawing

Parameter	Limit Values
Appearance(Ab097000)	Colourless liquid
IHC	Conform