

Shanghai Aladdin Biochemical Technology Co., LTD

Building No. 7, No. 11, Galileo Road, Pudong New Area, Shanghai 201210, China

Specification Sheet

Product Name:	Goat Anti-Chicken IgY H&L (HRP)
SKU #:	Ab137759
Version	3
Clone ID:	
Isotype:	IgG
Source:	Serum
Formulation:	Liquid
Concentration:	1 mg/mL
Buffer Solution:	0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide
Storage Conditions:	Store at -20°C

Parameter	Limit Values
Appearance(Ab137759)	Light brown to brown liquid
ELISA (1/20000)	≥1
Non-specific binding	Conform